

# Zoll X Series Transport Monitor/Defibrillator Checklist

Unit/Location \_\_\_\_\_ Device Number \_\_\_\_\_

Week of: \_\_\_/\_\_\_/\_\_\_ to \_\_\_/\_\_\_/\_\_\_

The defibrillator should be checked daily. Contact the **Service Center at 615.343.9600** to correct any issues found.

	Date	___/___/___	___/___/___	___/___/___	___/___/___	___/___/___	___/___/___
Initials							
1. Assure black circle <input checked="" type="checkbox"/> shows in <b>READY FOR USE</b> window.							
<ul style="list-style-type: none"> <li>If red slash circle <input type="checkbox"/> shows in the window, take the corrective action displayed (replace expired pads, connect to AC power, etc.)</li> <li>If after taking corrective action, black circle <input checked="" type="checkbox"/> does not show, take out of service and call the Service Center for replacement.</li> </ul>							
2. Assure defibrillator is connected to AC Power.							
3. Assure battery is fully charged as indicated when pressing the test button on the battery.							
4. Press the green power button and verify <b>SELF TEST PASSED</b> shows.							
<ul style="list-style-type: none"> <li>If <b>SELF TEST FAILED</b> shows, take defibrillator out of service and call the Service Center for replacement.</li> </ul>							
5. Disconnect defibrillator from A/C power.							
6. Connect therapy cable to test block.							
7. Press <b>ENERGY SELECT</b> buttons and set energy to 30 joules.							
8. Press the <b>CHARGE</b> button.							
9. When the Ready tone sounds, press <b>SHOCK</b> button until the shock is delivered.							
<ul style="list-style-type: none"> <li>If display shows <b>DEFIB SHORT TEST PASSED</b>, device is ready, push the green power button to power off.</li> <li>If display shows <b>DEFIB SHORT TEST FAILED</b>, take defibrillator out of service and call the Service Center for replacement.</li> </ul>							
10. Reconnect defibrillator to AC Power.							
11. Assure the following accessories are present:							
<ul style="list-style-type: none"> <li>Therapy cable</li> </ul>							
<ul style="list-style-type: none"> <li>Defibrillation Pads</li> </ul>							
<ul style="list-style-type: none"> <li>3-Lead Monitor Cable</li> </ul>							
<ul style="list-style-type: none"> <li>Monitoring Electrodes</li> </ul>							
<ul style="list-style-type: none"> <li>NIBP Cuff and Tubing</li> </ul>							
<ul style="list-style-type: none"> <li>Filterline H Set (ETCO2)</li> </ul>							

Signature	Initials
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____

Notes

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_