

Eliciting Epitope-Specific CD8⁺ T Cell Response by Immunization with Microbial Protein Antigens Formulated with α -Galactosylceramide: Theory, Practice, and Protocols

Pavlo Gilchuk, Frances C. Knight, John T. Wilson, and Sebastian Joyce

Abstract

CD8⁺ cytotoxic T lymphocytes confer protection against infectious diseases caused by viruses, bacteria, and parasites. Hence, significant efforts have been invested into devising ways to generate CD8⁺ T cell-targeted vaccines. Generation of microbe-free protein subunit vaccines requires a thorough knowledge of protective target antigens. Such antigens are proteolytically processed peptides presented by MHC class I molecules. To induce a robust antigen-specific CD8⁺ T cell response through vaccination, it is essential to formulate the antigen with an effective adjuvant. Here, we describe a versatile method for generating high-frequency antigen-specific CD8⁺ T cells through immunization of mice using the invariant natural killer T cell agonist α -galactosylceramide as the adjuvant.

Key words Adjuvant, α -Galactosylceramide, Antigen-specific CD8⁺ T cells, Microbial protein antigens, Mouse immunization

1 Introduction

Vaccination is one of medical science's greatest achievements; it has led to the eradication of several infectious diseases and to a reduction in morbidity and mortality from many others [1]. The majority of licensed vaccines consist of live, live-attenuated, or inactivated pathogens that were developed by empirical approaches [2]. Understanding vaccine-mediated correlates of immunity is critical for rational design of novel pathogen-free vaccines against infectious diseases for which no effective vaccine currently exists, or to improve the safety and/or efficacy of existing vaccines [2–6].

It is generally thought that the most effective vaccines confer protection via mediation of neutralizing antibodies, which are readily detectable by *in vitro* assays. Hence, most effort continues to be invested in designing vaccines that induce antibody-mediated

protective immunity [6, 7]. In recent years, however, the focus on developing T cell-targeted vaccines has increased, as it has become clear that successful vaccination relies on both antibody- and T cell-mediated immunity (reviewed in Ref. [3, 8]).

The indispensable role of cellular immunity was observed even prior to modern advances in microbiology and immunology, from the historically recorded cases known as “experiments of nature.” As early as the nineteenth century, Jenner and colleagues observed that a smallpox-immune mother had brought forth an infant who died a few days after birth owing to an infection contracted in utero, as the baby was born with pocks/pustules [9, 10]. Given that the mother’s antibody repertoire circulates within the fetus and babies are born with underdeveloped cellular immunity [11–13], we can now surmise that protection from smallpox requires cellular immunity. Another example of the critical role of cellular immunity includes a case report for progressive vaccinia virus (VACV) infection in a subject with normal anti-viral antibody responses [14].

Immune CD8⁺ T cells recognize infected cells via short microbial peptide epitopes presented by the host’s major histocompatibility complex (MHC)-encoded class I molecules. Of note, MHC class I molecules present cytoplasmic antigens because their assembly with peptides predominantly occurs within the endoplasmic reticulum, though under inflammatory conditions it can occur within vacuoles as well (reviewed in Refs. [15, 16]). Recent studies have shown that CD8⁺ T cells play a critical role in conferring protective immunity against a variety of infectious diseases, including those caused by respiratory viruses [17–21], poxviruses [22], cytomegalovirus [23], human immunodeficiency virus (HIV; [24]), *Plasmodium* species [25, 26], and *M. tuberculosis* [27]. Hence, T cell epitope-based vaccines can complement current efforts in developing safe and effective antibody-targeted subunit vaccines.

The challenge impeding the development of T cell-targeted subunit vaccines is identifying which of the microbe’s peptide epitopes confer protective immunity amongst the many that are presented [28]. Such epitopes can be discovered by indirect approaches that can ascertain protein subunit immunogenicity and protective capability in surrogate animal models [29, 30]. Accordingly, there is a need for the design and implementation of an effective vaccination strategy that will elicit sufficiently high-frequency, antigen-specific CD8⁺ T cells to allow thorough characterization of the response to subunit antigens. Unlike replicative microbial vaccines, heterologous proteins and peptides are poorly immunogenic when administered alone. Therefore, to induce a robust antigen-specific CD8⁺ T cell response without the use of infectious agents or antigen-loaded exogenous dendritic cells (DCs) [31–33], it is essential to formulate the antigen with an effective adjuvant [34]. Successful CD8⁺ T cell-targeted immunization has been demonstrated for peptide epitopes that are formulated with various toll-like receptor






(TLR) ligands [31, 35]. High-frequency CD8⁺ T cells can be induced by peptides formulated with polyinosinic-polycytidylic acid (poly(I:C); TLR3 agonist) and co-stimulatory anti-CD40 antibody [36]. Peptide immunogenicity can be improved via chemical engineering of peptides and adjuvants with amphiphilic albumin-binding tags [37]. The immunogenicity of peptides can also be modulated through engineered delivery systems; these are discussed elsewhere [38, 39] and, hence, not reviewed here.

It should be noted that full-length protein antigens are typically more desirable for vaccination studies than peptide epitopes. This is because immunization with proteins allows elicitation of both antibody- and T cell-mediated immune responses [28, 40–42]. In addition, one protein may contain several immune epitopes, which is a desirable feature for human vaccines [43]. Because MHC class I molecules present cytoplasmic antigens, the priming of naïve CD8⁺ T cells requires cross-presentation of exogenous antigens by antigen-presenting cells (APCs) that have taken up the administered protein. Small subsets of DCs—the tissue-resident, generally non-migratory, CD103⁺ as well as splenic CD8⁺ DCs in mice, or CD1c⁺ and CD141⁺ DCs in humans—play critical roles in cross-presentation of exogenous antigens (reviewed in Refs. [15, 16]). Antigen cross-presentation requires the activation and licensing of the appropriate DC subsets; CD4⁺ T cells are the classical helper cells for DC activation and licensing (reviewed in Ref. [44]). So also, CD4⁺ T cell help is essential for the differentiation of effector and memory CD8⁺ T cells (reviewed in Ref. [44]). Because MHC class II molecules—which also present processed peptide antigens generated in vacuolar compartments to control the functions of CD4⁺ T cells—are highly variable, vaccine design requires the identification of helper epitopes and their inclusion in CD8⁺ T cell-targeted vaccines. Hence, harnessing a non-polymorphic antigen-presenting molecule and invariant T cells as vaccine targets that can provide help can be a viable option for enhancing CD8⁺ T cell responses to subunit antigens. CD1d and the restricted semi-invariant natural killer T (NKT) cells have shown significant promise as experimental vaccine targets [44]; their biology has been recently reviewed elsewhere [45].

Immunogenicity of whole-protein antigens formulated with the NKT cell agonist α -galactosylceramide (α GalCer) as the adjuvant has demonstrated that this formulation stimulates robust antigen-specific T cell and antibody responses in mice [41, 46–48]. α GalCer is a synthetic glycolipid that closely resembles the marine sponge *Agelas mauritianus*-derived agelasphin 9b as well as self and bacterial glycolipid antigens (see Table 1 and references therein). Upon immunization, professional APCs, such as DCs, macrophages, and B cells, take up the administered protein antigen and α GalCer. Antigens undergo processing to short peptide epitopes that are presented by MHC class I molecules expressed





Table 1
Structure and properties of selected synthetic, microbial, and self NKT cell agonists

Lipid (class) ^a	Chain length ^b	Structure	Agonist ^c	References
αGalCer (GSL)	C18 C24:1		IFN-γ, IL-4 Self	[67]
Agel 9b (GSL)	C17 (C ₁₆ -Me) phyto C24		Anti-tumour; <i>Agelas mauritianus</i>	[68, 69]
αGalCer (GSL)	C18-phyto C26 <i>Common form used in functional studies (this is the molecule used in this chapter)</i>		Very strong Robust IFN-γ, IL-4, and other cytokines; synthetic analogue of Agel 9b (KRN7000)	[70]
DB06-1 (GSL)	C18-phyto C26 (C ₂ S)		Strong; strong IFN-γ; IL-10 upon re-activation	[71]
α-C-GalCer (GSL)	C18-phyto C26 1,1'-C-glycoside		Weak (mo ^d)-to-none (hu); IFN-γ; synthetic	[72]

EPH-αGalCer (GSL)	C18-phyto C10-fluoro-phenyl		Very strong IFN-γ; synthetic	[73]
NU-αGalCer (GSL)	C18-phyto C26 C6''-naphthylurea		Very strong Robust IFN-γ; synthetic	[74]
PyrC- αGalCer (GSL)	C18-phyto C26 C6''-(pyridin-4-yl)-carbamate		Very strong Robust IFN-γ, IL-12; synthetic	[75]
SMC124 (GSL)	C22-C ₁₁ -cyclo-propyl C26		Strong IFN-γ; synthetic	[76]
EF77 (GSL)	C18-phyto C21:1-C ₁₀ -cyclo-propyl		Strong IFN-γ; synthetic	[76]

(continued)

Table 1
(continued)

Lipid (class) ^a	Chain length ^b	Structure	Agonist ^c	References
OCH (GSL)	C9-phyto C24		Weak (mo)-to-none (hu) IL-4 (low-to-no IFN- γ); synthetic	[77]
C20-diene (GSL)	C18-phyto C20:2		Strong IL-4 (low-to-no IFN- γ); synthetic	[78]
α GalACer (GSL)	C18-phyto C14		Weak; <i>Sphingomonas</i> spp.	[79, 80]
Asp B (GSL)	C20:2-C ₉ Me C16-C ₂ OH		Weak <i>Aspergillus</i> <i>fumigatus</i>	[81]

^a *Agel* agelasphin, *Asp B* asparamide B, *GalCer* galactosylceramide, *GalACer* galacturonosylceramide, *GSL* glycosphingolipid, *IL* interleukin

^b Sphingosine/phytosphingosine chain length indicated first and N-acyl chain length second

^c Agonist strength based on Ref. [82]

^d Relative potency in comparison to α GalCer; *mo* mouse, *hu* human




Fig. 1 Strategy by which a robust antigen-specific CD8⁺ T cell response is elicited in mice immunized with antigenic protein formulated with the glycolipid adjuvant α GalCer. Upon immunization, professional antigen-presenting cells, such as dendritic cells (DCs), endocytose administered protein antigen and α GalCer. The antigen undergoes processing into short peptide epitopes that are presented by MHC class I molecules on DCs for recognition via the T cell receptor (TCR) of antigen-specific CD8⁺ T cells. α GalCer is presented by MHC class I-like CD1d molecules and recognized by the TCR of semi-invariant NKT cells. Recognition of α GalCer rapidly activates NKT cells, which produce cytokine/chemokine messengers, thereby trans-activating DCs and initiating a cross talk between members of the innate and adaptive immune systems. NKT cell licensing of DCs facilitates antigen cross presentation by a less well-understood mechanism(s) [42, 44, 65]. Together, this provides stimulus to naïve antigen-specific CD8⁺ T cell precursors that then proliferate and differentiate into effector and memory CD8⁺ T cells [45]

by APCs for recognition via the T cell receptor (TCR) of antigen-specific CD8⁺ T cells. The MHC class I-like, lipid antigen-presenting CD1d molecules present α GalCer and the TCR of NKT cells recognize the co-complex. Recognition of α GalCer rapidly activates NKT cells. Activated NKT cells produce cytokine/chemokine messengers, thereby trans-activating DCs and initiating a crosstalk between cells of the innate and adaptive immune systems. Activated NKT cells and DCs together stimulate naïve antigen-specific CD8⁺ T cell precursors, which then proliferate and differentiate into effector and memory CD8⁺ T cells (Fig. 1; [28, 44, 45, 49, 50]). By virtue of its ability to communicate to virtually all cell types of the immune system, α GalCer has proven to be a potent adjuvant that enhances antibody and CD4 T cell responses to foreign antigens (reviewed in Ref. [45]).

Several studies have reported that α GalCer is an effective adjuvant for the elicitation of epitope-specific CD8⁺ T cells against the model protein antigen ovalbumin [41, 51–53], or for enhancing antigen-specific CD8⁺ T cell responses upon immunization with live viral vaccines [49, 54]. We recently demonstrated the potency of α GalCer as a CD8⁺ T cell adjuvant for in vivo immunogenicity and protection studies which exploited pathogen-derived protein antigens ([28] and unpublished data). It is noteworthy that the activation of NKT cells by α GalCer results in the production of proinflammatory as well as immunoregulatory cytokines and chemokines (reviewed in [45]). Since DC activation and licensing as well as CD8⁺ T cell differentiation require proinflammatory cytokines and chemokines, efforts to synthesize and identify α GalCer analogues that elicit proinflammatory activity are being sought but are not discussed here (*see* Table 1; reviewed in Ref. [55]). Herein, we describe protocols for the preparation of protein antigens and α GalCer adjuvant, mouse immunization, and assessment of antigen-specific CD8⁺ T cell response, while analysis of antigen-specific antibody response and protective immunity is not described.

Experiments that assess immunogenicity require a large quantity of pure protein antigens, which, unlike model antigens, are not typically available commercially. Instead, the desired proteins are produced by recombinant DNA methods in *Escherichia coli* and then purified to a high degree of purity through engineered tags with the use of affinity resin. We designed, cloned, and produced nine recombinant VACV-derived protein subunits (rA3L_{62–319}, rD1R_{565–844}, rD5R_{330–470}, rE2L_{26–301}, rF4L_{1–319}, rJ6R_{188–466}, and rL4R_{33–249}, as well as epitope-engineered subunits L4R/B8R_{70–79} and L4R/A34R_{32–90}), all of which contain a C-terminal histidine-tag to facilitate purification (Fig. 2; [28]). Upon immunization of mice, α GalCer-formulated antigens induced robust epitope-specific CD8⁺ T cell and antibody responses that were detectable in blood and in various tissues (Fig. 3a–c). Vaccination with adjuvant and microbial antigen induced protective immunity against subsequent microbial challenge (Fig. 4 and unpublished data). Our immunization approach is suitable for elicitation of antigen-specific CD8⁺ T cell and antibody responses to recombinant proteins that are produced as insoluble inclusion bodies (IB). Results that have emerged from the use of approaches described herein demonstrate the efficacy of the proposed immunization strategy for generating high-frequency antigen-specific CD8⁺ T cells against microbial protein antigens that can be exploited for immunologic and vaccine studies.




Fig. 2 Design of recombinant VACV-derived protein antigen subunits that contain immune CD8⁺ T cell epitopes. Schematic chart of engineered proteins showing nine VACV antigens with amino acid positions selected for the design of a truncated recombinant subunit. Shown are the location and sequences of immune CD8⁺ T cell epitopes, and their N- and C-terminal tag sequences introduced to facilitate expression and purification of recombinant proteins in *E. coli*. Soluble protein rL4R₃₃₋₂₄₉ was engineered with B8R₇₀₋₇₉ and A34R₈₂₋₉₀ immune epitopes to facilitate expression of soluble antigens by *E. coli*. All CD8⁺ T cell epitopes are restricted to the HLA-B*07:02 (B7.2) molecules. This figure is adapted from our previous publication with permission from the Journal of Clinical Investigation [28]




Fig. 3 Efficient stimulation of antigen-specific immune responses is achieved through prime and boost immunization with α GalCer-formulated protein antigens. **(a)** Experimental design. Purified protein antigen is formulated with α GalCer and administered to mice by intraperitoneal or intranasal route as indicated. Antigen-specific immune responses were assessed ~6 days after booster immunization. **(b)** Representative example showing that immunization with α GalCer-formulated VACV-derived recombinant protein subunits elicited a robust epitope-specific CD8⁺ T cell response. Individual proteins containing HLA-B7.2 restricted CD8⁺ epitopes (see Fig. 2) were formulated with α GalCer and administered to groups of B7^{tg} mice IP as in shown in **(a)**. On d6 after boost, various tissues and blood were harvested from each group of immunized mice and assessed with B7.2




Fig. 4 Protection from lethal VACV infection upon vaccination with VACV-derived antigenic protein that contains CD8⁺ epitope. **(a)** Vaccination and challenge strategy. Three groups of mice ($n=5$ mice per group) were immunized by IP route with α GalCer-formulated VACV-derived protein antigens rJ6R₈₈₋₄₆₆ (contains protective CD8⁺ epitope), or engineered Δ rJ6R₃₀₃₋₃₁₁/L4R₃₇₋₄₅ (contains non-protective CD8⁺ epitope), or mock-vaccinated as indicated. Vaccinated mice were challenged with VACV and observed during the ensuing 18 days for severity of disease. Elicitation of J6R₃₀₃₋₃₁₁ (protective epitope) or L4R₃₇₋₄₅ (non-protective epitope) epitope-specific CD8⁺ was confirmed with corresponding tetramers in blood of vaccinated mice on d8 after second boost. **(b)** Protection of protein and α GalCer-vaccinated mice from lethal respiratory VACV challenge. All mock-vaccinated mice succumbed to the disease by d8 post-infection. Both protein and α GalCer-vaccinated groups were protected from severe weight loss and death. Replacing protective CD8⁺ epitope J6R₃₀₃₋₃₁₁ in rJ6R₈₈₋₄₆₆ with non-protective L4R₃₇₋₄₅ epitope showed that protection is mediated by J6R₃₀₃₋₃₁₁-specific CD8⁺ T cells with partial contribution from other immune responses against the protein antigen. *** $p < 0.001$; * $p < 0.05$ as compared to mock on d8 post challenge as defined by ANOVA with Dunnett's post-test; mean \pm SEM. Symbols indicate proteins used for vaccination

Fig. 3 (continued) tetramer for epitope-specific CD8⁺ T cells. Epitope-specific CD8⁺ T cells were detected using a dual-fluorochrome labeling approach as described previously [28, 66]. Contour plots are gated on live CD8⁺ T lymphocytes; numbers indicate percent epitope-specific CD8⁺ T cells. Panel was reproduced from Ref. 28 with permission from the Journal of Clinical Investigation. **(c)** Representative example showing that immunization with α GalCer-formulated VACV-derived recombinant protein subunits elicited antigen-specific antibody responses. SDS-PAGE (*left*) and western blot (*right*) of purified VACV-derived subunits. Each lane was loaded with $\sim 4 \mu\text{g}$ protein for SDS-PAGE and $\sim 0.5 \mu\text{g}$ protein for immunoblot. Blots were probed with 1:500 sera dilution from mice primed and boosted with the indicated individual proteins, and developed with anti-mouse HRP-conjugated secondary antibodies. *M* molecular weight standards, from top to bottom: 250, 150, 100, 75, 50, 37, 25, 20, 15 kDa

2 Materials

Reagents and equipment described herein can be substituted with equivalent alternatives from other vendors.

2.1 Adjuvant Preparation

1. Fisherbrand™ class A clear glass threaded 6 ml vials with caps (Fisher Scientific).
2. α GalCer (Funakoshi, Japan); *see* Table 1 for structure. To dissolve α GalCer, use 5.6 % sucrose, 0.75 % l-histidine, and 0.5 % Tween 20 in sterile water for injection (WFI) for cell culture with heating at 80 °C until dissolved (*see* Note 1). Make a 0.2 mg/ml stock solution and store at –20 °C.

2.2 Antigen Preparation

The immunization protocol reported herein was validated with the use of vaccinia virus (VACV)-derived recombinant proteins designed and produced in our laboratory [28] and with the model antigen ovalbumin. The following materials are required to obtain recombinant protein antigens for immunizing mice:

1. *E. coli* BL21-Gold(DE3) Competent Cells (Agilent Technologies) transformed with plasmid encoding desired protein antigen.
2. LB broth: To 800 ml H₂O add 10 g tryptone, 5 g yeast extract, and 10 g NaCl. Adjust pH to 7.5 with 0.5 M NaOH and bring the final volume to 1 L with dH₂O. Sterilize the solution by autoclaving.
3. Kanamycin (Km): Make 50 mg/ml stock solution in sterile dH₂O and store at –20 °C.
4. Glucose: Make 40 % (w/v) stock solution and sterilize by autoclaving.
5. LB agar plates containing 50 µg/ml Km and 1 % glucose.
6. Isopropyl β -d-thiogalactoside (IPTG): Make 1 mM stock in sterile dH₂O and store at –20 °C.
7. Sterile phosphate-buffered saline (PBS), 1×: 137 mM NaCl, 2.7 mM KCl, 10 mM Na₂HPO₄, 2 mM KH₂PO₄.
8. Ni-NTA agarose.
9. Cell lysis buffer A: 50 mM NaH₂PO₄, 300 mM NaCl, 5 mM imidazole (pH 8.0). Sterilize using Complete Filtration Unit with 0.2 µm filter (VWR) and store at room temperature (RT).
10. Lysozyme.
11. Wash buffer B: 50 mM NaH₂PO₄, 300 mM NaCl, 10 mM imidazole (pH 8.0). Filter sterilize (*see* Subheading 2.2, item 9) and store at room temperature (RT).

12. Elution buffer C: 50 mM NaH₂PO₄, 300 mM NaCl, 250 mM imidazole (pH 8.0). Filter sterilize (*see* Subheading 2.2, **item 9**) and store at RT.
13. Solubilization/wash buffer D: 20 mM Tris-HCl, 300 mM NaCl, 8 M urea, 10 mM imidazole (pH 8.0) (*see* **Note 2**). Filter sterilize (*see* Subheading 2.2, **item 9**) and store at RT.
14. Elution buffer E: 20 mM Tris-HCl, 300 mM NaCl, 8 M urea, 100 mM imidazole (pH 8.0). Filter sterilize (*see* Subheading 2.2, **item 9**) and store at RT.
15. 10–12 % Sodium dodecyl sulfate (SDS)-polyacrylamide gel (*see* **Note 3**).
16. Sodium dodecyl sulfate-polyacrylamide gel electrophoresis (SDS-PAGE) running buffer: 0.025 M Tris-HCl (pH 8.3), 0.192 M glycine, 0.1 % SDS.
17. Loading dye solution (5×): 0.03 M Tris-HCl (pH 6.8), 10 % SDS, 25 % 2-mercaptoethanol, 0.1 % bromophenol blue, 45 % glycerol. Store at –20 °C.
18. Bovine serum albumin (BSA) lyophilized powder: Prepare 0.2 mg/ml stock in 1× loading dye solution. Store at –20 °C.
19. Spectra Multicolor Broad Range Protein Ladder (Thermo Fisher).
20. GelCode Blue Stain Reagent (Thermo Fisher).
21. Optional: Pierce High Capacity Endotoxin Removal Spin Columns, 0.25 ml (Thermo Fisher).
22. Slide-A-Lyzer Dialysis Cassette, 10 K MWCO, 3 ml (Thermo Fisher).
23. Pierce disposable polypropylene columns, 5 ml (Thermo Fisher).
24. Falcon disposable polypropylene centrifuge tubes, 50 ml (Corning).
25. Disposable Safe-Lock Tubes, 1.5 ml (Eppendorf).
26. Disposable centrifuge tubes, 250 ml (Corning).
27. Nalgene baffled shake flask, 2000 ml (Sigma).
28. Disposable syringes, 5 ml (BD Biosciences).
29. Disposable Millex-GV Syringe Filter Unit, 0.22 μ m (EMD Millipore).
30. Disposable 28 mm syringe filter, 0.45 μ m (Corning).
31. Nalgene Oak Ridge High-Speed PPCO Centrifuge Tubes, 50 ml (Thermo Scientific).
32. Orbital shaker incubator.
33. 550 Sonic Dismembrator (Fisher Scientific).

34. Precision Microprocessor-Controlled 280 Series Water Bath (Thermo Scientific).
35. Mini-PROTEAN 3 Mini Vertical electrophoresis system (Bio-Rad) with Power Station 300 plus power supply (Labnet International, Inc).
36. Eppendorf 5415D microcentrifuge.
37. Refrigerated Sorvall RT7 centrifuge with RTH-750 rotor.
38. Sorvall RC-5B refrigerated centrifuge with SS-34 rotor.
39. DU 530 UV/Vis Spectrophotometer (Beckman).

2.3 Immunizations

1. Mice: All mice used for immunization were 6–10-week-old males from the following strains: B6-*KOD0;B*07;02tg* (B7^{tg}) human leukocyte antigen (HLA) class I transgenic (tg) mice were bred in house. B7^{tg} mouse is described elsewhere [56]. C57BL/6 (B6) mice were purchased from Jackson Laboratories (Bar Harbor, ME). All animal procedures were approved by the Institutional Animal Care and Use Committee (IACUC) at Vanderbilt University.
2. Anesthesia: 10 mg/ml KETAVED (ketamine HCl injection, USP; Vedco) and 1 mg/ml xylazine HCl Injection (xylazine; Vanderbilt pharmacy) in sterile WFI. Ketamine/xylazine anesthesia should be approved by an IACUC at user institution.
3. Purified recombinant antigen or albumin from chicken egg white (ovalbumin, Ova) at 0.5–2 mg/ml in sterile PBS (*see Note 4*).
4. Sterile insulin syringes, 1 ml (BD Biosciences).
5. Sterile filter pipet tips, 1–200 μ l (VWR).

2.4 Lymphocyte Isolation

1. SmartBox Auto CO₂ System for rodent euthanasia (Euthanex).
2. Refrigerated Sorvall RT7 centrifuge with RTH-750 rotor and plate adaptors.
3. Hausser Scientific Hemocytometer (Fisher Scientific).
4. Forceps.
5. Surgical scissors.
6. Falcon sterile disposable cell strainer, white, 70 μ m (Fisher Scientific).
7. BD sterile disposable Slip Tip syringes, 5 ml (Fisher Scientific).
8. Sterile disposable Falcon polypropylene centrifuge tubes, 15 ml (Fisher Scientific).
9. Falcon sterile disposable 60 mm \times 15 mm Petri dishes (Corning).
10. Optional: Capillary tubes, heparinized (Fisher Scientific) or plain (Fisher Scientific).

11. Dimethyl sulfoxide (DMSO).
12. Dasatinib (LC Laboratories): Dissolve to 1 mM in DMSO and store in 0.05 ml aliquots at -20°C .
13. FACS buffer: PBS containing 2% heat-inactivated fetal bovine serum (FBS) and 50 nM dasatinib. Make and use fresh (*see* **Note 5**).
14. Collagenase (Sigma).
15. Complete RPMI medium (cRPMI): RPMI-1640 with l-glutamine and HEPES supplemented with 10% FBS, $1\times$ penicillin/streptomycin solution and $1.75\ \mu\text{l}$ of cell culture grade 2-mercaptoethanol.
16. Lung digestion medium: cRPMI (*see* Subheading 2.4, **item 15**) containing 2 mg/ml collagenase and 50 nM dasatinib.
17. Gibco ACK lysing buffer (Thermo Fisher).
18. Trypan blue solution, 0.4% in PBS (Fisher Scientific).

2.5 Tetramer Staining, Intracellular Cytokine Staining (ICCS), and Flow Cytometry

1. 10 mM Peptide epitope (purity $\geq 95\%$; Schufer-N, Denmark) dissolved in sterile DMSO.
2. Brefeldin A (Sigma). Prepare stock solution at 10 mg/ml in methanol.
3. BD GolgiStop Protein Transport Inhibitor (BD Biosciences).
4. Phorbol 12-myristate 13-acetate (PMA) (InvivoGen): Make 1 mg/ml stock solution in sterile WFI and store at -20°C . Protect from light.
5. Ionomycin (Sigma): Dissolve in sterile DMSO to make 3 mM ionomycin solution and store at -20°C .
6. cRPMI (*see* Subheading 2.4, **item 15**).
7. FACS buffer (*see* Subheading 2.4, **item 13**).
8. BD Cytofix/Cytoperm Fixation and Permeabilization Solution (BD Biosciences).
9. BD Perm/wash buffer, $10\times$ solution (BD Biosciences).
10. Optional: AccuCheck counting beads (Thermo Fisher).
11. UltraComp eBeads (eBioscience).
12. Falcon 96 Well Round Bottom Cell Culture Plate (Corning).
13. Fisherbrand Polypropylene Microtiter Tubes, 1.2 ml (Fisher Scientific).
14. Refrigerated Sorvall RT7 centrifuge with RTH-750 rotor and plate adaptors.
15. LSR-II flow cytometer (BD Biosciences).
16. FlowJo flow cytometry software (Tree Star, Inc).
17. Tetramers, antibodies, and viability dyes: *See* Table 2.

Table 2
Flow cytometry reagents

Analysis ^a	Reagent	Antibody clone	Company ^c
Tet	B8R ₇₀₋₇₉ (FPKNDFVSF)/B7.2-PE	— ^b	—
Tet	L4R ₃₇₋₄₅ (FPRSMLSIF)/B7.2-PE	—	—
Tet	A3L ₁₉₂₋₂₀₀ (SPSNHHILL)/B7.2-PE	—	—
Tet	D1R ₆₈₆₋₆₉₄ (HPRHYATVM)/B7.2-PE	—	—
Tet	D1R ₈₀₈₋₈₁₇ (RPSTRNFFEL)/B7.2-PE	—	—
Tet	D5R ₃₇₅₋₃₈₃ (LPKEYSSEL)/B7.2-PE	—	—
Tet	J6R ₃₀₃₋₃₁₁ (MPAYIRNTL)/B7.2-PE	—	—
Tet	Ova ₂₅₇₋₂₆₄ (SIINFEKL)/H-2Kb-PE	—	—
Tet	Anti-CD8α-Pacific Blue	Clone 53-6.7	BD Biosciences
ICCS	Anti-CD8α-PerCP-Cy5.5	Clone 53-6.7	BD Biosciences
Tet	Anti-B220-FITC	Clone RA3-6B2	BD Biosciences
ICCS	Anti-B220-APC-Cy7	Clone RA3-6B2	BD Biosciences
Tet	Anti-CD11c-FITC	Clone N418	Tonbo Biosciences
ICCS	Anti-CD3e-PE	Clone 145-2C11	BD Biosciences
Tet	Anti-CD4-FITC	Clone H129.19	BD Biosciences
Tet	Anti-CD11b-FITC	Clone M1/70	Tonbo Biosciences
ICCS	Anti-CD107a-FITC	Clone 1D4B	BD Biosciences
ICCS	Anti-IFN-γ-allophycocyanin (APC)	Clone XMGI.2	BD Biosciences
Tet	Propidium iodide	—	BD Biosciences
ICCS	Fixable Viability Dye eFluor 450	—	eBioscience

^aTet p/MHC class I tetramer used for staining, ICCS intracellular cytokine staining

^bp/B7.2, H2K^b monomers and tetramers were made in-house according to published methods ([28, 83], see also “Production protocols” at <http://tetramer.yerkes.emory.edu/>); all tetramers are labeled with phycoerythrin (PE) and stored as 50 µg/ml stock at +4 °C

^cRepresentative vendors from which the reagents were purchased. These reagents are also available from the other vendors as well

3 Methods

Follow institutional biosafety and waste disposal regulations when working with chemicals and biologicals. Experimentation involving mice must comply with the experimenters’ IACUC regulations. The procedures described here were approved by IACUC at Vanderbilt University.

3.1 Preparation of Antigen

This section describes the preparation of recombinant protein antigen that was used to immunize mice. The sample protocol described below includes the production, purification, and formulation of a VACV protein produced as a recombinant protein by biosynthesis in *Escherichia coli* (*E. coli*). This protocol was validated for ten different recombinant proteins, three of which were produced in soluble form, with the remaining seven proteins produced in insoluble form as inclusion bodies (IB). Isolated protein antigens proved to be suitable for in vivo immunogenicity and protection studies. Each investigator should optimize this protocol to ensure optimal yield, purity, and solubility of the protein of interest.

1. Plate bacteria transformed with the desired plasmid on freshly prepared LB agar plate containing 50 $\mu\text{g}/\text{ml}$ Km and 1 % glucose; incubate overnight at +37 °C. Follow general rules of aseptic techniques when working with *E. coli*; that is, perform these procedures around the blue flame of a Bunsen burner.
2. Inoculate single colony into 2 ml of LB broth containing 50 $\mu\text{g}/\text{ml}$ Km and 1 % glucose; incubate overnight at +37 °C (*see Note 6*).
3. Inoculate 250 ml LB broth containing 50 $\mu\text{g}/\text{ml}$ Km (no glucose) with 1 ml of bacterial culture prepared as from Subheading 3.1, **step 2**, in 2000 ml Nalgene baffled shake flask; incubate at +37 °C with vigorous shaking (~150 rpm) to $A_{600\text{ nm}}=0.8\text{--}1.0$ (this level of growth takes ~4 h).
4. Induce expression of the recombinant gene by adding IPTG to 1 mM and continue incubation for an additional 4 h. Optional: Perform analysis of cellular fractions, i.e., total cell lysate as well as soluble and insoluble fractions, by 10–12 % SDS-PAGE to verify expression and solubility of the recombinant protein. Pause point: Bacterial culture can be stored overnight at +4 °C, or if the protein is produced as IB, cells can be harvested by centrifugation as in Subheading 3.1, **step 5**, and stored frozen at –80 °C, if necessary, for a long period of time.
5. Pellet cells by centrifugation at $3300\times g$ for 20 min at +4 °C. Discard the supernatant; carefully resuspend the pellet in 5 ml of cell lysis buffer A; then add 5 ml of the same buffer containing 4 mg/ml lysozyme; transfer to 50 ml polypropylene tube and incubate for 30 min on ice.
6. Sonicate with the 550 Sonic Dismembrator ten times for 10 s on ice using a setting of 5 on the Amplitude Control Knob, which is the highest setting for the micro-tip probe. Follow the manufacturer's instructions for the instrument in use.
7. Transfer the cell lysate into 50 ml Nalgene Oak Ridge High-Speed PPCO Centrifuge tube, and centrifuge at $7670\times g$ for 20 min at +4 °C. Harvest the pellet and proceed to **step 8** if protein is produced as IB; harvest supernatant and proceed to **step 9** if soluble protein is produced.

8. To purify IB, discard supernatant and resuspend the pellet in 10 ml cell lysis buffer A without lysozyme by sonication as in Subheading 3.1, **step 6**, followed by centrifugation as in **step 7**. To solubilize protein in the IB, resuspend pellet in 0.5 ml PBS by sonication as in Subheading 3.1, **step 6** first, then add 9.5 ml solubilization/wash buffer D (*see Note 7*).
9. Filter the supernatant through a 45 μm filter into a 15 ml polypropylene tube; add 2 ml Ni-NTA slurry, mix, and incubate with rotation for 2 h at +4 °C (*see Note 8*).
10. Pellet the resin at 300 $\times g$ for 5 min at +4 °C and carefully remove the supernatant (*see Note 8*). Resuspend resin in 14 ml of the appropriate wash buffer: for native protein purification, use wash buffer B; for solubilized denatured protein, use solubilization/wash buffer D. Centrifuge the solution at 300 $\times g$ for 5 min at +4 °C. Repeat this wash step.
11. Re-suspend the resin in 5 ml of the corresponding wash buffer as in Subheading 3.1, **step 10** and transfer to a disposable 5 ml polypropylene column. Wash the column four times with 5 ml wash buffer B (for native protein purification) or D (for solubilized denatured protein purification).
12. Elute with 3 ml elution buffer C (for native protein purification) or elution buffer E (for denatured protein purification). Sterilize the eluted protein by filtration through a Millex 0.22 μm filter and store at +4 °C on ice. Use for immunization within the next 2 days.
13. Assess purity of isolated protein using 10–12% SDS-PAGE. SDS-PAGE analysis has been described elsewhere [57] and, hence, is not described here. Quantify purified protein; one approach is to use a known amount of BSA protein that is separated by SDS-PAGE on the same gel as the recombinant proteins for quantification by densitometry. An average yield of purified recombinant VACV protein is ~1.5–6 mg, which corresponds to ~0.5–2 mg/ml in the eluted fraction.
14. Recommended: Follow vendor's protocol to further purify protein with the use of endotoxin removal spin column (Subheading 2.2, **item 21**) or, alternatively, using the recently published protocol (*see Note 9*). Otherwise, proceed to **step 15** (optional) or **step 16**.
15. Optional: Dialyze protein that was eluted under native conditions to exchange the elution buffer with PBS (*see Note 10*).
16. Remove urea from protein that was purified under denaturing conditions (*see Note 11*). Dialyze eluted protein against 500 ml PBS at RT; this will promote protein precipitation when urea is removed (*see Note 12*). Centrifuge at 10,000 $\times g$ for 5 min at RT, remove supernatant, and re-solubilize precipi-

tate in 300 μ l of 1% SDS. This should yield ~5–20 mg/ml of protein. Incubate at +37 °C until completely dissolved. Store at RT and use for immunization over the next 2 days. Note that SDS precipitates upon thawing SDS-solubilized protein stored at –20 °C, as will the protein stored with it.

3.2 Immunizing Mice

The most common routes for immunizing a mouse are intraperitoneal (IP), subcutaneous (SQ), intravenous (IV), and intranasal (IN). SQ is commonly used for tracking CD8⁺ T cell responses in the draining lymph node [37]. Intramuscular (IM) administration is not generally recommended, as the muscle of the mouse is small. The route of immunization used determines trafficking properties of antigen-specific CD8⁺ T cells and influences the anamnestic immune response at the initial site of pathogen invasion [58]. Hence, choice of immunization route will depend on the investigator's goal and infection model. Below, we describe protocols for IP and IN immunizations, which have been shown to elicit robust and protective CD8⁺ T cell responses against lethal VACV infection in mice ([28] and unpublished data). IP immunization induces a robust systemic response, resulting in antigen-specific CD8⁺ T cells circulating through the blood and secondary lymphoid organs, such as the spleen. IN immunization, in addition to eliciting a systemic response will stimulate a robust local mucosal immune response within the lungs [59]. Generating a robust antigen-specific CD8⁺ T cell response generally requires two immunizations—prime and boost. Experimentation involving mice must comply with the investigator's IACUC regulations.

1. Thaw 0.2 mg/ml α GalCer stock solution at RT, and mix gently by pipetting (*see* **Note 13**).
2. Prepare appropriate volume of inoculum by mixing protein with α GalCer. For IN immunization, prepare 50 μ l of mixture that contains 20–100 μ g of protein (45 μ l of 0.5–2 mg/ml solution) and 1 μ g of α GalCer (5 μ l of 0.2 mg/ml stock) per mouse. For IP immunization, prepare 200 μ l of mixture that contains the same amount of protein and α GalCer per mouse as for IN immunization. If protein is insoluble and contains SDS, adjust SDS concentration with sterile PBS accordingly to final SDS concentration \leq 0.05% in the inoculum. For example, add 4 μ l re-solubilized protein at 5–20 mg/ml in 1% SDS from Subheading 3.1, **step 16**, and 5 μ l of 0.2 mg/ml α GalCer stock to 191 μ l of PBS (*see* **Note 14**). Proceed with **step 3** for IP immunization, or **step 4** for IN immunization.
3. For IP immunization, inject 200 μ l of the inoculum IP into a B7.2^{tg} mouse (*see* **Note 15**).
4. For IN immunization, anesthetize a B7.2^{tg} mouse with 150 μ l of ketamine-xylazine (*see* **Note 16**). After the mouse is deeply anesthetized (confirmed by the absence of foot reflexes in

response to a mild pinch), gradually release 50 μ l of the inoculum into both nostrils at the same time with a micropipette. Adjust the rate of release so as to allow the mouse to inhale the inoculum without forming bubbles. This should cause a rapid increase in the breathing rate. Hold the mouse in an upside down position for another couple of minutes until its breathing gradually returns to normal. Observe mice until recovery from anesthesia (*see* **Note 17**).

5. Boost mice 14 days after primary immunization as follows: prepare mixture of protein and α GalCer as in **step 1** and inoculate mice as in **step 2** (*see* **Notes 18–20**).

3.3 Isolation of Lymphocytes and Flow Cytometry

Note that all experimentation involving mice must comply with the investigator's IACUC regulations. In addition, flow cytometry should be performed by trained personnel and in accordance with institutional regulations. Refer elsewhere for detailed overviews of flow cytometry acquisition and data analysis [**60, 61**].

1. Sacrifice mouse using SmartBox Auto CO₂ System.
2. For each mouse, prepare two 60 mm Petri dishes and place on ice. For IP-immunized mice, remove the spleen and put on top of a 70 μ m cell strainer placed inside the Petri dish. For IN-immunized mice, remove the spleen and lungs and place inside separate Petri dishes, each with a 70 μ m cell strainer. Harvest spleen or lungs from non-immunized mice as negative controls. Proceed with **step 3** to isolate splenic lymphocytes or **step 6** to isolate lung lymphocytes (*see* **Note 21**).
3. To each spleen, add 2.5 ml ACK lysis buffer at RT. Using the plunger end of a 5 ml syringe, mash the tissue over the cell strainer to disperse leukocytes. Pass the cell suspension one more time through the cell strainer and transfer into a 15 ml polypropylene centrifuge tube. Incubate for 2 min at RT to lyse red blood cells (RBCs). Add 10 ml cRPMI with dasatinib to dilute ACK lysis buffer.
4. Centrifuge the cell suspension at 300 $\times g$ for 5 min at +4 °C. Decant supernatant.
5. Resuspend each splenocyte suspension by gentle pipetting in 5 ml FACS buffer and place on ice. To stain splenocytes with various fluorochrome-conjugated antibodies, proceed to **step 9** (*see* **Note 22**).
6. Rinse the harvested lungs with PBS to remove blood. Mince tissue with a scalpel; transfer the minced lung tissue into a 15 ml polypropylene centrifuge tube containing 2 ml cRPMI supplemented with dasatinib and 2 mg/ml collagenase and let digest for 1 h at +37 °C.

7. Pour the digested lung slices onto a 70 μ m cell strainer placed inside a Petri dish. Using the plunger end of a 5 ml syringe, mash the tissue over the cell strainer to disperse leukocytes.
8. Follow Subheading 3.3, steps 3–4, to lyse RBCs, prepare single-cell suspension, then resuspend cells from each lung in 2 ml of FACS buffer, and place on ice (*see Note 22*).
9. Count viable cells with a hemocytometer or counting beads. The expected recovery is $\sim 0.5\text{--}1 \times 10^7$ leukocytes per 1 ml of cell suspension from processed lungs or spleen (*see Note 23*).
10. Add 200 μ l of cells from each sample, including control non-immunized mouse leukocytes, into a 96-well round-bottom plate. Optional: If absolute count of tetramer-positive CD8⁺ T cells is required, we recommend using counting beads (*see Note 23*). Centrifuge plate at $828 \times g$ for 2 min at +4 °C. Discard the supernatant.
11. Resuspend the cells in 100 μ l of FACS buffer containing anti-B220-fluorescein isothiocyanate (FITC), -CD4-FITC, -CD11b-FITC, -CD11c-FITC (all at 1:200 dilution), -CD8 α -Pacific Blue (1:100 dilution), and 1.5 μ g/ml of corresponding R-phycoerythrin (PE)-labeled tetramer. Incubate for 1 to 2 h at +4 °C (*see Note 24*).
12. Dilute the staining reaction with 200 μ l of FACS buffer and spin plate at $800 \times g$ for 2 min at +4 °C to wash the cells. Discard the supernatant.
13. Resuspend the cells in 200–300 μ l of FACS buffer containing 1:2000 dilution of 50 μ g/ml propidium iodide stock solution, which incorporates into the nuclei of dead cells and, thereby, allowing differentiation between dead and viable cells in flow cytometry experiments. Transfer cells into 1.2 ml Microtiter Tubes; store on ice for up to 2 h before acquisition.
14. Prepare an acquisition template using FACSDiva software containing four fluorochromes to reflect the staining panel.
15. Set up compensation by using an unstained control and single stained compensation control samples.
16. Set up an acquisition layout containing the following plots: forward scatter (FSC)-A versus side scatter (SSC)-A. Create a gate for the lymphocyte population, and a FSC-W versus FSC-H followed by SSC-W versus SSC-H gate to discriminate singlet population. Create a histogram to display live/dead cells for propidium iodide with a gate to discriminate propidium iodide-negative live population; FITC (“dump” channel) versus CD8 α -Pacific Blue to create CD8⁺ T cell gate; CD8 α versus tetramer to create tetramer⁺ gate (Fig. 5).
17. Collect 5000–10,000 total CD8 α ⁺ events and export FCS files for the analysis in FlowJo software.




Fig. 5 Analysis of epitope-specific CD8⁺ T cells by tetramer staining. **(a)** Representative example showing gating strategy and identification of epitope-specific CD8⁺ T cells. Cells harvested from the lungs of rL4R/B8R₇₀₋₇₉ protein/αGalCer primed and boosted B7¹⁹ mice (IN route) were analyzed by four-color flow cytometry following staining with B8R₇₀₋₇₉/B7.2 tetramer. **(b)** Representative flow cytometry contour plots showing specificity of B8R₇₀₋₇₉/B7.2 tetramer staining upon immunization with αGalCer-formulated protein antigen

18. In FlowJo, to determine the frequency of antigen-specific, i.e., tetramer⁺, CD8⁺ T cells, apply gates as in Subheading 3.3, **step 16**. To determine background from non-specific tetramer staining, create tetramer gate for “CD8 α versus tetramer” subset for cells from non-immunized control mouse, then apply it to the cells from a mouse immunized with protein and adjuvant. To determine the frequency of tetramer⁺ CD8⁺ T cells, subtract the background from each experimental value. If using counting beads, refer to vendor’s protocol to determine absolute number for the desired lymphocyte subset.
19. Optional: The remaining unstained cells can be cryopreserved at -80°C in cRPMI containing 10% DMSO for CD8⁺ T cell phenotyping in additional experiments if needed at a later time (*see* **Note 25**).

3.4 CD8 T Cell Re-stimulation and Staining for Intracellular Interferon- γ

Staining for intracellular interferon- γ (IFN- γ) can be optionally used as a complementary assay to validate the presence of functional antigen-specific CD8⁺ T cells after immunization.

1. Harvest spleen and isolate lymphocytes as reviewed in Subheading 3.3, **steps 1–4**. Work under aseptic conditions using a laminar flow hood; do not use dasatinib in any medium during tissue processing and cell incubation. Include dasatinib in FACS buffer only during the staining procedure for flow cytometry.
2. Resuspend cells from each spleen in 5 ml of cRPMI. Add 200 μl of cell suspension containing $\sim 2 \times 10^6$ lymphocytes and 50 μl of anti-CD107a-FITC antibody (optional, *see* **Note 26**) to each of three wells of 96-well round-bottom plate. To the first well, add 50 μl of cRPMI containing 50 μM antigenic peptide epitope to 10 μM final. Add 50 μl of cRPMI to the second well; this will serve as a negative control for CD8⁺ T cell re-stimulation. Add 50 μl of cRPMI containing PMA and ionomycin to a final concentration of 50 ng/ml PMA plus 2 $\mu\text{g}/\text{ml}$ ionomycin to the third well; this will serve as a positive control for CD8⁺ T cell stimulation.
3. Incubate for 2 h at $+37^{\circ}\text{C}$ in CO_2 incubator (5% CO_2), then add 50 μl cRPMI containing a combination of vesicular transport inhibitors, such as brefeldin A (10 $\mu\text{g}/\text{ml}$ final) and GolgiStop (1:1500); incubate for an additional 4 h, for a total of 6 h. These inhibitors will allow for intracellular accumulation of secreted proteins including cytokines such as IFN- γ at levels that are readily detectable by flow cytometry.
4. Spin plate at $828 \times g$ for 2 min at $+4^{\circ}\text{C}$. Discard the supernatant, resuspend in 250 μl of PBS, and spin again at $828 \times g$ for 2 min at $+4^{\circ}\text{C}$ to wash the cells. Discard the supernatant.

5. Resuspend the cells in 100 μ l of PBS containing eFluor 450 amino-reactive viability dye (1:2000). Incubate 30 min at +4 °C.
6. Add 200 μ l of FACS buffer and centrifuge plate at 828 $\times g$ for 2 min at +4 °C to wash the cells. Discard the supernatant.
7. Resuspend the cells in 100 μ l of FACS buffer supplemented with dasatinib and containing anti-B220-APC-Cy7 (allophycocyanin-cyanine7 tandem dye) (1:100), -CD3 ϵ -PE (1:100) and -CD8 α -PerCP-Cy5.5 (peridinin chlorophyll-cyanine5.5 tandem dye) (1:200) antibodies (*see* **Note 27**). Incubate for 1 h at +4 °C.
8. Add 200 μ l of FACS buffer and spin plate at 828 $\times g$ for 2 min at +4 °C to wash the cells. Discard the supernatant.
9. Resuspend the cells in 100 μ l of BD Cytofix/Cytoperm and incubate for 10 min at +4 °C.
10. Wash the cells twice by adding 200 μ l of 1 \times BD Perm/wash and centrifuging at 828 $\times g$ for 2 min at +4 °C. Discard the supernatant.
11. Resuspend the cells in 100 μ l of 1 \times BD Perm/wash containing anti-IFN- γ -allophycocyanin (1:200) antibody. Incubate for 1 h at +4 °C.
12. Add 200 μ l of 1 \times BD Perm/wash and centrifuge plate at 828 $\times g$ for 2 min at +4 °C to wash the cells. Discard the supernatant, resuspend in 200–300 μ l of FACS buffer, and transfer cells into 1.5 ml microtubes. Stained cells can be stored on ice for up to 4 h before data acquisition. Alternatively, cells can be fixed as in Subheading 3.4, **step 9**, washed with FACS buffer as in **step 8**, resuspended in 200–300 μ l of FACS buffer, and stored at +4 °C up to 16 h before data acquisition. Longer storage times have not been tested.
13. Prepare an acquisition template for six fluorochromes and set up compensation as detailed in Subheading 3.3, **steps 14–15**.
14. Set up an acquisition layout containing the following plots: forward scatter (FSC)-A versus side scatter (SSC)-A. Create a gate for the lymphocyte population and an SSC-W versus SSC-H gate to identify the singlet population. Create a histogram to display live/dead cells for eFluor450 with a gate to identify live cells; B220-APC-Cy7 versus CD3 ϵ -PE to create T cell gate; CD8 α -PerCP-Cy5.5 versus CD3 ϵ -PE to create CD8 $^{+}$ T cell gate; CD107a-FITC versus IFN- γ -allophycocyanin to create CD107a $^{+}$ IFN- γ $^{+}$ CD8 $^{+}$ T cell gate (Fig. 6).
15. Analyze sample as detailed in Subheading 3.3, **steps 17–18**. To determine frequency of CD107a $^{+}$ IFN- γ $^{+}$ CD8 $^{+}$ T cells, subtract the background from a well without antigenic peptide epitope from each experimental value.




Fig. 6 Analysis of epitope-specific CD8⁺ T cells by intracellular cytokine staining (ICCS). **(a)** Representative example showing gating strategy and identification of epitope-specific CD8⁺ T cells. Cells harvested from the spleen of protein/ α GalCer primed and boosted mice (IP route) were re-stimulated *in vitro* with the indicated antigenic peptide epitope, then stained with the indicated antibodies and analyzed by six-color flow cytometry. **(b)** Representative flow cytometry plots showing that IFN- γ was specifically induced by a subset of CD8⁺ T cells after peptide epitope re-stimulation

4 Notes

1. For preparation and storage of α GalCer, use only glass vials to minimize leaching of impurities out of a plastic container and adsorption of the lipid to plastic (*see* Subheading 2.1, **item 1**).
2. All urea-containing buffers should be made fresh. At physiological pH, urea produces cyanate that can cause carbamylation of proteins by reacting with amino, carboxyl, and sulfhydryl groups, which may alter their stability and function.

3. Refer to methods described elsewhere [57] for making protein electrophoresis gel and performing SDS-PAGE.
4. Proteins purified from IBs may contain up to 0.05 % SDS for solubility.
5. Low concentrations of dasatinib were reported to substantially improve CD8⁺ T cell detection with tetramers by preventing T cell receptor internalization [62]. Dasatinib is a reversible inhibitor that targets a wide variety of protein kinase C family enzymes [63]. Hence, we routinely use dasatinib in FACS buffer instead of sodium azide to maintain the TCR at the cell surface and for improved staining and cell viability.
6. Use freshly transformed bacterial colonies grown on plates; or alternatively, plate transformed bacteria from frozen glycerol stock.
7. Prior resuspension of IB in PBS through vigorous sonication facilitates further solubilization of proteins within IB by using urea-containing solubilization/wash buffer D. Centrifuge cells at $7670 \times g$ for 20 min at +4 °C and harvest the supernatant.
8. Save an aliquot of the supernatant for SDS-PAGE analysis.
9. LPS is a TLR4 ligand that could potentiate T cell and antibody responses and, hence, we recommend its removal. Nonetheless, we did not find trace amounts of endotoxin contaminants and other impurities co-purifying with the protein preparations to contribute substantially to the elicitation of a CD8⁺ T cell response in the presence of α GalCer [28]. In addition, LPS removal in the past has resulted in dramatic losses of VACV protein. Should it become necessary to remove co-purifying endotoxins, use the recently described endotoxin clearance method by including a non-ionic detergent wash step during metal chelation affinity purification of the protein antigen [64].
10. This step is optional because components of elution buffer C do not affect the performance of the immunization protocol or antigen-specific immune responses.
11. Proteins purified under denaturing conditions contain a high concentration of urea, which is incompatible with the immunization protocol. Removal of urea by common laboratory approaches may result in protein aggregation. Stimulation of antigen-specific T cell responses *in vivo* does not require native conformation of protein antigen, but does require the protein to be in a reasonably soluble form. Although soluble urea-free protein can be obtained through *in vitro* refolding, such a procedure is laborious and requires optimization for each protein. We found that a low SDS concentration (0.01–0.05 %) is compatible with IP immunization protocol. This residual low concentration of SDS can promote protein solubility after the removal of urea by dialysis.

12. Some proteins can be refolded successfully and remain soluble after dialysis.
13. α GalCer solution should be clear; if appears cloudy, re-heat at +80 °C until dissolved. Repeated freeze and thaw does not affect the reagent's performance as long as the solution upon thawing remains clear and has not changed color.
14. Prepare mixture of protein with α GalCer right before immunization; usage of SDS-containing protein for IN immunization has not been validated.
15. For the first experiment, we recommend setting up a control wherein the mouse is inoculated IP with the same formulation but without protein antigen. Leukocytes harvested from this mouse after immunization will serve as a control to determine the background of the tetramer staining assay. An additional control for the performance of the reagents and immunization protocol that can be considered is the inclusion of a C57BL/6 mouse immunized by either the IP or IN route with 1 μ g of α GalCer and 50 μ g of ovalbumin (Ova) as detailed in the main protocol. Ova-specific CD8⁺ should be readily detected with the readily available H2K^b/SIINFEKL tetramer (NIH Tetramer Core, Emory University, Atlanta, USA; <http://tetramer.yerkes.emory.edu/>) in the spleen on day 6–14 after boost.
16. Depending on the weight of the mouse, the amount of anesthesia can vary from 125 to 165 μ l.
17. For the first experiment, we recommend setting up a control wherein the mouse is inoculated IN with the same formulation containing PBS and α GalCer but without protein antigen. Lymphocytes harvested from this mouse after immunization will serve as a control to determine the background of the tetramer staining assay.
18. Usually, antigen-specific CD8⁺ T cells can be detected at low frequency on days 8–14 after primary immunization. At least one booster immunization substantially increases their frequency, which is more desirable for immunological and vaccination studies.
19. For the first time, before the mouse is sacrificed and leukocytes harvested, it may be necessary to confirm that immunization has induced antigen-specific CD8⁺ T cells of detectable frequency. This can be accomplished by assessing mouse blood leukocytes. A robust antigen-specific CD8⁺ T cell response should be readily detected systemically the blood on day 6–14 after boost.
20. Although it is not reviewed in the protocol here, protein antigen and α GalCer immunization also induces antigen-specific antibody responses that can be assessed in mouse serum (but *see* Fig. 3C).

21. Harvested organs can be left on ice while additional mice are being processed or reagents are being prepared for up to 24 h. Longer storage times have not been tested.
22. Cell suspensions can be left on ice for up to 8 h while additional tissues are being processed. Longer storage times have not been tested.
23. We found that precise counting of CD8⁺ T cells in the lungs requires counting beads. Counting beads are added directly to staining solution to account for cell losses, assuming proportional bead and cell losses during washes. When working with counting beads, follow vendor's protocol at https://tools.thermofisher.com/content/sfs/manuals/PCB100_accu-check_beads_man.pdf.
24. The investigator should determine optimal staining conditions, such as concentration of tetramer, incubation temperature, and time; each of these parameters will differentially influence experimental outcome.
25. Cryopreservation does not affect staining performance but may affect viability. An expected lymphocyte viability after thawing is 80% or higher.
26. CD107a is a marker for degranulation activity of immune effector cells including CD8⁺ T cells. Staining with anti-CD107a antibody is optional and performed during re-stimulation of cells.
27. Upon in vitro re-stimulation, epitope-specific CD8⁺ T cells markedly down-regulate their T cell receptor from the cell surface and, hence, stain poorly with the complex of peptide/MHC tetramers.

Acknowledgements

Supported by Vanderbilt University Discovery Grant as well as VA Merit Award (BX001444) and NIH Contracts (AI040079), Research (AI042284, HL121139), Core (CA068485, DK058404), and Center (CA068485) grants.

References

1. Plotkin SA, Plotkin SL (2011) The development of vaccines: how the past led to the future. *Nat Rev Microbiol* 9:889–893
2. Zepp F (2010) Principles of vaccine design—lessons from nature. *Vaccine* 28(Suppl 3):C14–C24
3. Amanna IJ, Slifka MK (2009) Wanted, dead or alive: new viral vaccines. *Antiviral Res* 84:119–130
4. Bevan MJ (2011) Understand memory, design better vaccines. *Nat Immunol* 12:463–465
5. Pulendran B, Ahmed R (2011) Immunological mechanisms of vaccination. *Nat Immunol* 13:509–517
6. Koff WC, Gust ID, Plotkin SA (2014) Toward a human vaccines project. *Nat Immunol* 15:589–592
7. Amanna IJ, Slifka MK (2011) Contributions of humoral and cellular immunity to vaccine-induced protection in humans. *Virology* 411:206–215

8. Moss B (2011) Smallpox vaccines: targets of protective immunity. *Immunol Rev* 239:8–26
9. Rees AR (2015) The antibody molecule: from antitoxins to therapeutic antibodies, 1st edn. Oxford Medical Histories/Oxford University Press, Oxford, p 384
10. Jenner E (1809) Two cases of Small-Pox Infection communicated to the Foetus in Utero under peculiar circumstances, with additional remarks. *Med Chir Trans* 1:271–277
11. Adkins B, Leclerc C, Marshall-Clarke S (2004) Neonatal adaptive immunity comes of age. *Nat Rev Immunol* 4:553–564
12. Swamy GK, Wheeler SM (2014) Neonatal pertussis, cocooning and maternal immunization. *Expert Rev Vaccines* 13:1107–1114
13. Verhasselt V (2015) Is infant immunization by breastfeeding possible? *Philos Trans R Soc Lond B Biol Sci*. doi:10.1098/rstb.2014.0139
14. O'Connell CJ, Karzon DT, Barron AL, Plaut ME, Ali VM (1964) Progressive vaccinia with normal antibodies. A case possibly due to deficient cellular immunity. *Ann Intern Med* 60:282–289
15. Blum JS, Wearsch PA, Cresswell P (2013) Pathways of antigen processing. *Annu Rev Immunol* 31:443–473
16. Segura E, Amigorena S (2015) Cross-presentation in mouse and human dendritic cells. *Adv Immunol* 127:1–31
17. Thomas PG, Keating R, Hulse-Post DJ, Doherty PC (2006) Cell-mediated protection in influenza infection. *Emerg Infect Dis* 12:48–54
18. Brown LE, Kelso A (2009) Prospects for an influenza vaccine that induces cross-protective cytotoxic T lymphocytes. *Immunol Cell Biol* 87:300–308
19. Kohlmeier JE, Woodland DL (2009) Immunity to respiratory viruses. *Annu Rev Immunol* 27:61–82
20. Kremer M, Suezer Y, Volz A, Frenz T, Majzoub M, Hanschmann KM, Lehmann MH, Kalinke U, Sutter G (2012) Critical role of perforin-dependent CD8⁺ T cell immunity for rapid protective vaccination in a murine model for human smallpox. *PLoS Pathog* 8:e1002557
21. Goulding J, Bogue R, Tahiliani V, Croft M, Salek-Ardakani S (2012) CD8⁺ T cells are essential for recovery from a respiratory vaccinia virus infection. *J Immunol* 189:2432–2440
22. Gordon SN, Cecchinato V, Andresen V, Heraud JM, Hryniewicz A, Parks RW, Venzon D, Chung HK, Karpova T, McNally J, Silvera P, Reimann KA, Matsui H, Kanehara T, Shimura Y, Yokote H, Franchini G (2011) Smallpox vaccine safety is dependent on T cells and not B cells. *J Infect Dis* 203:1043–1053
23. Li CR, Greenberg PD, Gilbert MJ, Goodrich JM, Riddell SR (1994) Recovery of Hla-restricted cytomegalovirus (Cmv)-specific T-cell responses after allogeneic bone-marrow transplant - correlation with Cmv disease and effect of ganciclovir prophylaxis. *Blood* 83:1971–1979
24. Goulder PJR, Watkins DI (2008) Impact of MHC class I diversity on immune control of immunodeficiency virus replication. *Nat Rev Immunol* 8:619–630
25. Schmidt NW, Butler NS, Badovinac VP, Harty JT (2010) Extreme CD8⁺ T cell requirements for anti-malarial liver-stage immunity following immunization with radiation attenuated sporozoites. *PLoS Pathog* 6:e1000998
26. Epstein JE, Tewari K, Lyke KE, Sim BKL, Billingsley PF, Laurens MB, Gunasekera A, Chakravarty S, James ER, Sedegah M, Richman A, Velmurugan S, Reyes S, Li M, Tucker K, Ahumada A, Ruben AJ, Li T, Stafford R, Eappen AG, Tamminga C, Bennett JW, Ockenhouse CF, Murphy JR, Komisar J, Thomas N, Loyevsky M, Birkett A, Plowe CV, Loucq C et al (2011) Live attenuated malaria vaccine designed to protect through hepatic CD8(+) T cell immunity. *Science* 334:475–480
27. Woodworth JSM, Behar SM (2006) Mycobacterium tuberculosis-specific CD8(+) T cells and their role in immunity. *Crit Rev Immunol* 26:317–352
28. Gilchuk P, Spencer CT, Conant SB, Hill T, Gray JJ, Niu X, Zheng M, Erickson JJ, Boyd KL, McAfee KJ, Oseroff C, Hadrup SR, Bennink JR, Hildebrand W, Edwards KM, Crowe JE Jr, Williams JV, Buus S, Sette A, Schumacher TN, Link AJ, Joyce S (2013) Discovering naturally processed antigenic determinants that confer protective T cell immunity. *J Clin Invest* 123:1976–1987
29. Gerdts V, Littel-van den Hurk SV, Griebel PJ, Babiuk LA (2007) Use of animal models in the development of human vaccines. *Future Microbiol* 2:667–675
30. Gerdts V, Wilson HL, Meurens F, van Drunen Littel-van den Hurk S, Wilson D, Walker S, Wheler C, Townsend H, Potter AA (2015) Large animal models for vaccine development and testing. *ILAR J* 56:53–62
31. van Hall T, van der Burg SH (2012) Mechanisms of peptide vaccination in mouse models: tolerance, immunity, and hyperreactivity. *Adv Immunol* 114:51–76
32. Koup RA, Douek DC (2011) Vaccine design for CD8⁺ T lymphocyte responses. *Cold Spring Harb Perspect Med* 1:a007252
33. Remakus S, Rubio D, Ma X, Sette A, Sigal LJ (2012) Memory CD8⁺ T cells specific for a

- single immunodominant or subdominant determinant induced by peptide-dendritic cell immunization protect from an acute lethal viral disease. *J Virol* 86:9748–9759
34. Coffman RL, Sher A, Seder RA (2010) Vaccine adjuvants: putting innate immunity to work. *Immunity* 33:492–503
 35. Duthie MS, Windish HP, Fox CB, Reed SG (2011) Use of defined TLR ligands as adjuvants within human vaccines. *Immunol Rev* 239:178–196
 36. Cho HI, Celis E (2009) Optimized peptide vaccines eliciting extensive CD8⁺ T-cell responses with therapeutic antitumor effects. *Cancer Res* 69:9012–9019
 37. Liu H, Moynihan KD, Zheng Y, Szeto GL, Li AV, Huang B, Van Egeren DS, Park C, Irvine DJ (2014) Structure-based programming of lymph-node targeting in molecular vaccines. *Nature* 507:519–522
 38. Saroja C, Lakshmi P, Bhaskaran S (2011) Recent trends in vaccine delivery systems: a review. *Int J Pharm Invest* 1:64–74
 39. Moon JJ, Huang B, Irvine DJ (2012) Engineering nano- and microparticles to tune immunity. *Adv Mater* 24:3724–3746
 40. Hermans IF, Silk JD, Gileadi U, Salio M, Mathew B, Ritter G, Schmidt R, Harris AL, Old L, Cerundolo V (2003) NKT cells enhance CD4⁺ and CD8⁺ T cell responses to soluble antigen in vivo through direct interaction with dendritic cells. *J Immunol* 171:5140–5147
 41. Singh N, Hong S, Scherer DC, Serizawa I, Burdin N, Kronenberg M, Koezuka Y, Van Kaer L (1999) Activation of NK T cells by CD1d and alpha-galactosylceramide directs conventional T cells to the acquisition of a Th2 phenotype. *J Immunol (Cutting Edge)* 163:2373–2377
 42. Semmling V, Lukacs-Kornek V, Thaïss CA, Quast T, Hochheiser K, Panzer U, Rossjohn J, Perlmutter P, Cao J, Godfrey DI, Savage PB, Knolle PA, Kolanus W, Forster I, Kurts C (2010) Alternative cross-priming through CCL17-CCR4-mediated attraction of CTLs toward NKT cell-licensed DCs. *Nat Immunol* 11:313–320
 43. Oseroff C, Peters B, Pasquetto V, Moutaftsi M, Sidney J, Panchanathan V, Tschärke DC, Maillere B, Grey H, Sette A (2008) Dissociation between epitope hierarchy and immunoprevalence in CD8⁺ responses to vaccinia virus western reserve. *J Immunol* 180:7193–7202
 44. Gottschalk C, Mettke E, Kurts C (2015) The role of invariant natural killer T cells in dendritic cell licensing, cross-priming, and memory CD8(+) T cell generation. *Front Immunol* 6:379
 45. Hill TM, Bezbradica JS, Van Kaer L, Joyce S (2016) CD1d-restricted natural killer T cells. *Encyclopaedia Live Sci*. doi:10.1002/9780470015902.a0020180.pub2
 46. Fujii S, Shimizu K, Smith C, Bonifaz L, Steinman RM (2003) Activation of natural killer T cells by alpha-galactosylceramide rapidly induces the full maturation of dendritic cells in vivo and thereby acts as an adjuvant for combined CD4 and CD8 T cell immunity to a coadministered protein. *J Exp Med* 198:267–279
 47. Galli G, Pittoni P, Tonti E, Malzone C, Uematsu Y, Tortoli M, Maione D, Volpini G, Finco O, Nuti S, Tavarini S, Dellabona P, Rappuoli R, Casorati G, Abrignani S (2007) Invariant NKT cells sustain specific B cell responses and memory. *Proc Natl Acad Sci U S A* 104:3984–3989
 48. Kamijuku H, Nagata Y, Jiang X, Ichinohe T, Tashiro T, Mori K, Taniguchi M, Hase K, Ohno H, Shimaoka T, Yonehara S, Odagiri T, Tashiro M, Sata T, Hasegawa H, Seino KI (2008) Mechanism of NKT cell activation by intranasal coadministration of alpha-galactosylceramide, which can induce cross-protection against influenza viruses. *Mucosal Immunol* 1:208–218
 49. Reilly EC, Thompson EA, Aspeslagh S, Wands JR, Elewaut D, Brossay L (2012) Activated iNKT cells promote memory CD8⁺ T cell differentiation during viral infection. *PLoS One* 7, e37991
 50. Carreno LJ, Kharkwal SS, Porcelli SA (2014) Optimizing NKT cell ligands as vaccine adjuvants. *Immunotherapy* 6:309–320
 51. Ko SY, Ko HJ, Chang WS, Park SH, Kweon MN, Kang CY (2005) alpha-galactosylceramide can act as a nasal vaccine adjuvant inducing protective immune responses against viral infection and tumor. *J Immunol* 175:3309–3317
 52. Courtney AN, Thapa P, Singh S, Wishahy AM, Zhou D, Sastry J (2011) Intranasal but not intravenous delivery of the adjuvant alpha-galactosylceramide permits repeated stimulation of natural killer T cells in the lung. *Eur J Immunol* 41:3312–3322
 53. Nishimura T, Kitamura H, Iwakabe K, Yahata T, Ohta A, Sato M, Takeda K, Okumura K, Van Kaer L, Kawano T, Taniguchi M, Nakui M, Sekimoto M, Koda T (2000) The interface between innate and acquired immunity: glycolipid antigen presentation by CD1d-expressing dendritic cells to NKT cells induces the differentiation of antigen-specific cytotoxic T lymphocytes. *Int Immunol* 12:987–994

54. Guillonneau C, Mintern JD, Hubert FX, Hurt AC, Besra GS, Porcelli S, Barr IG, Doherty PC, Godfrey DI, Turner SJ (2009) Combined NKT cell activation and influenza virus vaccination boosts memory CTL generation and protective immunity. *Proc Natl Acad Sci U S A* 106:3330–3335
55. Duwaerts CC, Gregory SH (2011) Targeting the diverse immunological functions expressed by hepatic NKT cells. *Expert Opin Ther Targets* 15:973–988
56. Alexander J, Oseroff C, Sidney J, Sette A (2003) Derivation of HLA-B*0702 transgenic mice: functional CTL repertoire and recognition of human B*0702-restricted CTL epitopes. *Hum Immunol* 64:211–223
57. Green M, Sambrook J. (2012) Molecular cloning. A laboratory manual (Fourth Edition) (Cold Spring Harbor Laboratory Press), Vol 3, p Protocol 8
58. Belyakov IM, Ahlers JD (2009) What role does the route of immunization play in the generation of protective immunity against mucosal pathogens? *J Immunol* 183:6883–6892
59. Li AV, Moon JJ, Abraham W, Suh H, Elkhader J, Seidman MA, Yen M, Im EJ, Foley MH, Barouch DH, & Irvine DJ (2013) Generation of effector memory T cell-based mucosal and systemic immunity with pulmonary nanoparticle vaccination. *Sci Transl Med* 5: 204ra130
60. Sharrow SO (2002) Overview of flow cytometry. *Curr Protoc Immunol* Chapter 5: Unit 5 1
61. Roederer M (2002) Multiparameter FACS analysis. *Curr Protoc Immunol* Chapter 5: Unit 5 8
62. Lissina A, Ladell K, Skowera A, Clement M, Edwards E, Seggewiss R, van den Berg HA, Gostick E, Gallagher K, Jones E, Melenhorst JJ, Godkin AJ, Peakman M, Price DA, Sewell AK, Wooldridge L (2009) Protein kinase inhibitors substantially improve the physical detection of T-cells with peptide-MHC tetramers. *J Immunol Methods* 340:11–24
63. Montero JC, Seoane S, Ocana A, Pandiella A (2011) Inhibition of SRC family kinases and receptor tyrosine kinases by dasatinib: possible combinations in solid tumors. *Clin Cancer Res* 17:5546–5552
64. Zimmerman T, Petit Frere C, Satzger M, Raba M, Weisbach M, Dohn K, Popp A, Donzeau M (2006) Simultaneous metal chelate affinity purification and endotoxin clearance of recombinant antibody fragments. *J Immunol Methods* 314:67–73
65. Wagner CS, Grotzke JE, Cresswell P (2012) Intracellular events regulating cross-presentation. *Front Immunol* 3:138
66. Hadrup SR, Bakker AH, Shu CJ, Andersen RS, van Veluw J, Hombrink P, Castermans E, Thor Straten P, Blank C, Haanen JB, Heemskerk MH, Schumacher TN (2009) Parallel detection of antigen-specific T-cell responses by multidimensional encoding of MHC multimers. *Nat Methods* 6:520–526
67. Kain L, Webb B, Anderson BL, Deng S, Holt M, Costanzo A, Zhao M, Self K, Teyton A, Everett C, Kronenberg M, Zajonc DM, Bendelac A, Savage PB, Teyton L (2014) The identification of the endogenous ligands of natural killer T cells reveals the presence of mammalian alpha-linked glycosylceramides. *Immunity* 41:543–554
68. Morita M, Motoki K, Akimoto K, Natori T, Sakai T, Sawa E, Yamaji K, Koezuka Y, Kobayashi E, Fukushima H (1995) Structure-activity relationship of alpha-galactosylceramides against B16-bearing mice. *J Med Chem* 38:2176–2187
69. Natori T, Koezuka Y, Higa T (1993) Agelasphins, Novel α -galactosylceramides from the marine sponge *Agelas mauritianus*. *Tetrahedron Lett* 34:5591–5592
70. Kawano T, Cui J, Koezuka Y, Toura I, Kaneko Y, Motoki K, Ueno H, Nakagawa R, Sato H, Kondo E, Koseki H, Taniguchi M (1997) CD1d-restricted and TCR-mediated activation of α 14 NKT cells by glycosylceramides. *Science* 278:1626–1629
71. Birkholz AM, Girardi E, Wingender G, Khurana A, Wang J, Zhao M, Zahner S, Illarionov PA, Ueno X, Li M, Yuan W, Porcelli SA, Besra GS, Zajonc DM, Kronenberg M (2015) A novel glycolipid antigen for NKT cells that preferentially induces IFN- γ production. *J Immunol* 195:924–933
72. Schmieg J, Yang G, Franck RW, Tsuji M (2003) Superior protection against malaria and melanoma metastases by a C-glycoside analogue of the natural killer T cell ligand α -Galactosylceramide. *J Exp Med* 198:1631–1641
73. Li X, Fujio M, Imamura M, Wu D, Vasan S, Wong CH, Ho DD, Tsuji M (2010) Design of a potent CD1d-binding NKT cell ligand as a vaccine adjuvant. *Proc Natl Acad Sci U S A* 107:13010–13015
74. Aspeslagh S, Li Y, Yu ED, Pauwels N, Trappeniers M, Girardi E, Decruy T, Van Beneden K, Venken K, Drennan M, Leybaert L, Wang J, Franck RW, Van Calenbergh S, Zajonc DM, Elewaut D (2011) Galactose-modified iNKT cell agonists stabilized by an induced fit of CD1d prevent tumour metastasis. *EMBO J* 30:2294–2305

75. Aspeslagh S, Nemcovic M, Pauwels N, Venken K, Wang J, Van Calenbergh S, Zajonc DM, Elewaut D (2013) Enhanced TCR footprint by a novel glycolipid increases NKT-dependent tumor protection. *J Immunol* 191:2916–2925
76. Tyznik AJ, Farber E, Girardi E, Birkholz A, Li Y, Chitale S, So R, Arora P, Khurana A, Wang J, Porcelli SA, Zajonc DM, Kronenberg M, Howell AR (2011) Glycolipids that elicit IFN- γ -biased responses from natural killer T cells. *Chem Biol* 18:1620–1630
77. Miyamoto K, Miyake S, Yamamura T (2001) A synthetic glycolipid prevents autoimmune encephalomyelitis by inducing TH2 bias of natural killer T cells. *Nature* 413:531–534
78. Yu KO, Im JS, Molano A, Dutronc Y, Illarionov PA, Forestier C, Fujiwara N, Arias I, Miyake S, Yamamura T, Chang YT, Besra GS, Porcelli SA (2005) Modulation of CD1d-restricted NKT cell responses by using N-acyl variants of α -galactosylceramides. *Proc Natl Acad Sci U S A* 102:3383–3388
79. Kinjo Y, Wu D, Kim G, Xing GW, Poles MA, Ho DD, Tsuji M, Kawahara K, Wong CH, Kronenberg M (2005) Recognition of bacterial glycosphingolipids by natural killer T cells. *Nature* 434:520–525
80. Mattner J, Debord KL, Ismail N, Goff RD, Cantu C 3rd, Zhou D, Saint-Mezard P, Wang V, Gao Y, Yin N, Hoebe K, Schneewind O, Walker D, Beutler B, Teyton L, Savage PB, Bendelac A (2005) Exogenous and endogenous glycolipid antigens activate NKT cells during microbial infections. *Nature* 434:525–529
81. Albacker LA, Chaudhary V, Chang YJ, Kim HY, Chuang YT, Pichavant M, DeKruyff RH, Savage PB, Umetsu DT (2013) Invariant natural killer T cells recognize a fungal glycosphingolipid that can induce airway hyperreactivity. *Nat Med* 19:1297–1304
82. Joyce S, Girardi E, Zajonc DM (2011) NKT cell ligand recognition logic: molecular basis for a synaptic duet and transmission of inflammatory effectors. *J Immunol* 187:1081–1089
83. Rodenko B, Toebes M, Hadrup SR, van Esch WJ, Molenaar AM, Schumacher TN, Ovaas H (2006) Generation of peptide-MHC class I complexes through UV-mediated ligand exchange. *Nat Protoc* 1:1120–1132