

Vanderbilt - Fixed Medical Device Dressing Selection Guide

Objective: Protect tissue and minimize friction, shear and moisture from fixed devices using Mepilex® and Mepilex® Lite



Mepilex® Lite

PMM
4x4 #75441



Mepilex®


PMM
4x4 #75440
4x8 #83549

Notations

- ✦ Wear time: Up to 7 days
- ✦ Fenestrate/cut product PRN to accommodate tube sites
- ✦ When cutting products, leave backing film in place. Cut to desired shape
- ✦ Products listed on this guide are not suitable for fixation of life sustaining devices
- ✦ Dressings with Safetac® technology **DO NOT** require use of skin barrier products



CPAP/BiPAP



O2 Mask Strap




ET Tube with Tape




Nasal Cannulae



Brace/Cast



Tracheostomy Tie



Tracheostomy




Nasal Cannulae with
Ear Protection



C-Collar



Multipodis Boot






G-Tube

See reverse side for Mepilex®
Cutting Guide

Cutting Guide

Select the dressing size appropriate to cover the affected area and cut to customize around unique patient needs/body contours, device shapes.



The suggested topical management options and change rates are the treatment choice of Vanderbilt Medical Center and may not reflect the opinions of Mölnlycke Health Care.

The information provided herein is not to be construed as the practice of medicine or substituted for the independent medical judgment of a patient's treating clinician. This information, including but not limited to suggestions for product wear time, product selection and suggested use is based on generalizations and does not consider the unique characteristics of an individual's wound. Each patient's clinician shall remain solely responsible for assessing the severity of patient wounds, determining the appropriate treatment, and managing treatment of the wound. For additional information, please refer to the applicable product insert or contact Mölnlycke Health Care at 1-800-843-8497.

The Mölnlycke Health Care, Mepilex and Safetac, names and respective logos are registered trademarks of Mölnlycke Health Care AB.
Distributed by Mölnlycke Health Care US, LLC, Norcross, Georgia 30092. © 2014. Mölnlycke Health Care AB. 1-800-843-8497

MHC-2014-2088 Fac Rev 04.14